



# Disinfection Procedure

Document method code: DP 03 toilet

Date: 06/05/2020

Revision: 2

issue no: 2

Review dates: 10-2021

Page 1 of 1



## Disinfection of toilet

<b>equipment/area to be disinfected:</b>		toilet
<b>area:</b>		ablutions
<b>assigned to:</b>		disinfecting team/cleaning personnel
<b>chemical to be used:</b>		vitrodx®-4
<b>dilution (use ratio):</b>		<b>viruses:</b> 50ml per litre (1:20) <b>microbiological:</b> 20ml per litre (1:50)
requirements		procedure
<b>safety: mild irritant eyes:</b> wash with cold running water for 5 minutes <b>prolonged skin contact:</b> wash area with cold running water inhalation: remove to fresh air	<b>ppe:</b> goggles, mask, nitrile gloves	<b>surface disinfection</b> a) apply solution by low-pressure sprayer/ toilet brush, disposable cloth and allow to be rinsed with the next flush  
<b>rinse water temperature:</b>	ambient	
<b>cleaning equipment:</b>	toilet brush, low-pressure sprayer, disposable cloth or sponge	
<b>rinsing equipment:</b>	clean disposable cloth or sponge	
<b>drying equipment:</b>	natural ventilation	
<b>frequency: daily</b> (directly proportional to traffic)	allow surface contact for minimum of 10 mins (no flush)	
key inspection points:		
1. under rim 2. water surface		
<b>inspected by:</b>	<b>supervisor:</b>	