



Disinfection Procedure

Document method code: DP 08 floor

Date: 06/05/2020

Revision: 2

issue no: 2

Review dates: 10-2021

Page 1 of 1



Disinfection of floors

equipment/area to be disinfected:		floor
area:		general
assigned to:		disinfecting team/cleaning personnel
chemical to be used:		b bioactive™-4
dilution (use ratio):		viruses: 50ml per litre (1:20) microbiological: 20ml per litre (1:50)
requirements		procedure
safety: mild irritant eyes: wash with cold running water for 5 minutes prolonged skin contact: wash area with cold running water inhalation: remove to fresh air	ppe: goggles, mask, nitrile gloves	<p>Disinfection:</p> <ul style="list-style-type: none"> a) clean floor with a suitable detergent b) apply solution with a bucket and mop and allow recommended contact time c) rinse with clean water <p>Maintenance:</p> <ul style="list-style-type: none"> a) apply with low-pressure sprayer, evenly directly in front of mop with sweeping motion
rinse water temperature:	ambient	
cleaning equipment:	low-pressure sprayer, bucket and mop	
rinse equipment:	clean water and mop	
contact time:	5-15 minutes	
drying equipment:	natural ventilation	
frequency: daily [directly proportional to foot traffic in the area]	viruses: hourly for high traffic, daily for low traffic microbiological: every 4-5 hours	
key inspection points:		
1. floor surface (edges)		
inspected by:	supervisor:	

